Brass Credit Program





2A Warehouse and Capital Cartridge are happy to offer our customers additional savings by providing a **Brass Credit Program** that can be applied to future orders on our website. It is extremely simple to do and will save money on your next purchase.

- 1) Create an account on 2AWarehouse.com or CapitalCartridge.com, depending on where you want the credit applied.
- 2) Web store credit only provided for spent cases made out of brass. No calibler restrictions.
- 3) No steel or aluminum cases accepted. Make sure to remove shotgun hulls, loose primers, and any live rounds. If you send these in, the weight will be subtracted.
- 4) Sorting is NOT needed (mixed calibers are ok our machines will do the sorting)
- 5) Bag brass inside of shipper box
- 6) Include a completed Customer Credit form in each shipper box

NO restictions on amount of brass you can send in.

- 7) Send in USPS Large Flate Rate Box (Use extra tape especially at corners and seams)
- 8) Brass credit will be applied to your account. Credit can be used on anything on our website including new ammunition.
- 9)The credit must be applied to your account prior to placing your order. When you place your order, you must mark the check box to apply store credit.

Example: Your order is \$1,000.00, send us \$999.00 in brass using our Brass Credit Program. After you apply your store credit, you will only be charged \$1.00 plus shipping! It's that easy!

Once we have received your brass, a store credit will be applied to your account. Credits typically take 1-3 business days to post after receipt of your package.

	ut at line, place bottom portion in each box of brass you are shipping)
Weight:	Number of Boxes:
Name:	
Phone Number:	
Email Address:	
	ne email used to create your account on 2AWarehouse.com or CapitalCartridge.com)
Total Number of Box	es you will be sending:
Total Weight:	
What Website To Ap	ply Credit To (2A Warehouse or Capital Cartridge):

Shipping Address:

Capital Cartridge- BC 9702 Wallisville Rd Suite B **Admin Use Only**

Total Weight Measured: